

Genesol 37

Strongly Acidic Membrane Cleaner

Genesol 37 has been developed as a general-purpose cleaner for inorganic scale and iron deposits found in most types of membrane system. It is also highly effective at removing organic deposits from severely fouled systems.

Application

Genesol 37 is effective when used as a 1-2% solution and alternately soaked and recycled through the membrane system for 1-2 hours.

Cleaning solution pH should be the minimum recommended by the membrane manufacturer. If the pH starts to rise during the cleaning process additional Genesol 37 should be added or the cleaning solution replaced.

After cleaning, Genesol 37 cleaning solutions should be slowly neutralised to pH 6-7 using an alkaline solution before discharge, as the presence of acid solutions can be damaging to concrete and cement lined discharge pipes. Before returning the plant to service, the membrane system should be thoroughly flushed, until the pH of the feed and the concentrate are similar. The product water conductivity may initially be high but will return to normal after a few hours of operation.

Genesol 37 is fully miscible with water in all proportions.

Health and Safety

Genesol 37 is a mixture of nitric and phosphoric acids and care should be taken during handling. Solutions of Genesol 37 are corrosive to carbon steel and many other metals and materials of construction, so all chemical tanks, pumps and pipe-work must be made from acid-resistant materials.


A material safety data sheet is available on request.

Classification

CHIP: Corrosive | UN: Corrosive

Packaging

Available in 25 kg kegs, 250 kg drums and 1,300 kg IBCs. The shelf life is two years under normal storage conditions.



Removes Inorganic Scale



Removes Iron Deposits

- ✓ Compatible with Reverse Osmosis (RO) and Nanofiltration (NF) membranes
- ✓ Effective against inorganic scale and iron deposits
- ✓ Effective against proteinous deposits
- ✓ Wide application of use including the food process industry
- ✓ Liquid product, easy to use
- ✓ Low dose rate ensures cost effective treatment

Typical properties

Appearance:	pale liquid
pH as supplied:	< 1
Specific gravity:	1.30 - 1.34
Freezing point:	< - 5°C

The information provided in this data sheet is believed to be true and accurate. Genesys International Ltd. accepts no product liability as the use of its products are outside the company's control.